

 **tekno valves**
50 YEARS OF EXCELLENCE

Post-Type & Toggle Operated Pin Index Valves for Medical Gases

Detailed Series Catalogue



BPTN-12 & BPKN-12



Your safety is valued

ISO 9001 & TPED certified valve manufacturer



	Page
1. Operating Principle & Identifying Features	2
1. Features & Benefits for Best-in-Class Performance	
Post-Type Valve	3.1
Toggle Operated Valve	3.2
3. Material of Construction & Assembly Arrangement	
Post-Type Valve	4.1
Toggle Operated Valve	4.2
4. Disassembly, Inspection & Assembly Instructions	5
5. Product Selection Guide – Valve Item Code Matrix	6
5. List of Approved Gases	7



Series BPKN-12 & BPTN-12

Identifying Features

BPKN-12 is key (post) & BPTN-12 is toggle operated pin index soft seated valves with two-piece stem construction. The design is used for medical gas cylinders with Yoke connection for filling & discharge.

The upper & lower stem interface with a square drive. The threads are located on the lower stem & the upper stem is free-floating. The design uses O-ring to create a seal around the upper stem. PEEK thrust washer is capsuled with the upper stem & acts as an anti-friction ring when the upper stem rotates to open & close the valve under high pressure. Any leakage past the packing nut metallic sealing due to external impact on the packing nut is prevented by packing nut O-ring. The lower stem has PA 66 seat insert to ensure low torque closure. The design does not use any fluorinated polymer whose by-products of combustion are toxic.

Valves are MRI approved up to 3 Tesla as per ASTM F2052-15. Markings, where provided, to indicate safety in the Magnetic Resonance Environment is as per ASTM F2503-13. *

Recommended Opening Procedure

It is recommended that the valves always be opened gradually in anticlockwise direction. Opening the valve fully causes the lower stem to ride upwards on its threads until it contacts the upper stem. Valves in fully open position can be mistaken as closed by inexperienced or untrained operators. The operator should always check the valve's position by attempting to close the valve, never by trying to open the valve.

Recommended Closing Procedure

Close the cylinder valve by rotating the key / toggle in the clockwise direction.

Valve Installation

1. Follow ISO 13341, where applicable, for valving procedure & torque guidelines.
2. For NGT threads, use hand tight + 3 turns wrench tight to install the valve in the cylinders
(Refer https://drive.google.com/file/d/1E0H1B_Z4rBb7ddQJ6R897duZPmFSzHCH/view?usp=sharing)
3. For UNF threads, use valve installation torque as below:

Inlet connection	Nm	ft.lb
0.750-16 UNF-2A	80-100	60-75
1.125-12 UNF-2A	100-130	75-95

CAUTION

1. Ensure valving tools (e.g. sockets or jaws) are not in contact with the PRD.
2. Valves should not be over-torqued into the cylinder to avoid high stresses leading to overload failures, especially in parallel thread valves.
3. Do not overtorque the valve in open or closed direction.
4. Do not attempt to replace soft seat in the lower stem.
5. Repair & maintenance shall be carried out by trained personnel.
6. Yokes used for filling & discharge shall only come in contact with the intended sealing surface.


* Markings indicate safety of the valves in Magnetic Resonance Environment and it should not be construed to be applicable for the full cylinder package. It is the responsibility of the user to ensure that the cylinder as well as other fittings are also approved and marked for MRI environment.



Features & Benefits for Best-in-Class Performance

Series BPTN-12 (Toggle Operated Valve)

Valve shown with Taper Inlet & CG-1 PRD



All dimensions are in mm

Dimensions shown are for 17E inlet

- a Depends upon inlet connection
- b Also available in short toggle (length – 34.5)

Design Specifications		
	Metric	English
Minimum life	2000 cycles	
Pressure rating	250 bar	3625 psig
Oxygen pressure surge test	20 cycles at 250 bar	20 cycles at 3625 psig
Temperature range	-46 °C to +85 °C	-51 °F to +185 °F
Flow coefficient (C _v)	0.16	
Minimum closing torque	2 Nm	1.5 ft.lb
Packing nut installation torque	50 Nm	37 ft.lb
PRD installation torque*		
- CG-1	16 Nm	12 ft.lb
- CG-4	7 Nm	5 ft.lb
Lubricant	Gleitmo 599	
MAX weight of cylinder package mass for which valve can be used without protection	22.2 kg	49 lb


*Optional

Testing & Certification

- Valves meet EN ISO 10297:2017 & CGA V-9:2019
- Valves are certified by BAM to European Transportable Pressure Equipment Directive (TPED) & available with Π mark
- All valves are MRI approved up to 3 Tesla as per ASTM F2052-15 & stamped as per ASTM F2503-13
- PRD complies with CGA S-1.1
- Production testing as per EN ISO 14246/CGA V-9



Series BPKN-12 (Post-type Valve)
Valve shown with Parallel Inlet & CG-4 PRD



All dimensions are in inch

Dimensions shown are for 3/4-16 UNF 2A inlet

a Depends upon inlet connection

Design Specifications		
	Metric	English
Minimum life	2000 cycles	
Pressure rating	250 bar	3625 psig
Oxygen pressure surge test	20 cycles at 250 bar	20 cycles at 3625 psig
Temperature range	-46 °C to +85 °C	-51 °F to +185 °F
Flow coefficient (C _v)	0.16	
Minimum closing torque	2 Nm	1.5 ft.lb
Packing nut installation torque	50 Nm	37 ft.lb
PRD installation torque*		
- CG-1	16 Nm	12 ft.lb
- CG-4	7 Nm	5 ft.lb
Spindle failure torque	22-24 Nm	16-18 ft.lb
Lubricant	Gleitmo 599	
MAX weight of cylinder package mass for which valve can be used without protection	22.2 kg	49 lb

*Optional



Testing & Certification

- Valves meet EN ISO 10297:2017 & CGA V-9:2019
- Valves are certified by BAM to European Transportable Pressure Equipment Directive (TPED) & available with **T** mark
- All valves are MRI approved up to 3 Tesla as per ASTM F2052-15 & stamped as per ASTM F2503-13
- PRD complies with CGA S-1.1
- Production testing as per EN ISO 14246/CGA V-9



Material of Construction & Assembly Arrangement


Series BPKN-12 (Post-type Valve) Valve shown with Parallel inlet and CG-4 PRD



Material of Construction & Assembly Arrangement



Series BPTN-12 (Toggle Operated Valve) Valve shown with Taper inlet and CG-1 PRD





Disassembly, Inspection & Assembly Instructions

Series BPTN-12 & BPKN-12

Disassembly of Valve

1. Place the valve in a vice with soft jaws. The holding fixture must securely grip the valve body (1) on the wrench flats so that there is no damage to the valve body surface, internal bores, PRD, inlet and outlet threads.
2. Using a 16 mm socket wrench, unscrew the packing nut (4e/5) in counter clockwise direction.
For post-type valve, upper stem assembly (4), i.e., upper stem (4d), O-ring (4b), Back-Up Ring (4c) and thrust washer (4a) will remove with the packing nut (5).
For toggle operated valves, packing nut assembly (4) toggle (4g) with pin (4f), upper stem, O-ring, Back-Up Ring and thrust washer will remove with the packing nut (4e).
3. Use the upper stem square to remove the lower stem assembly (3) from the valve chamber, by rotating it counter clockwise.
4. Remove the PRD (5/9) (if necessary) using 9.50 mm socket wrench or HEX box wrench (for CG-1) or T-30 TORX star tool (for CG-4) by rotating in counter clockwise direction.

Inspection of Valve & Components

1. Valve body (1)
 - a. Blow, using dry compressed air or Nitrogen, valve body (1) internal surfaces to remove any dirt or foreign particles.
 - b. Inspect the valve body for seat damage and thread wear.
 - c. Inspect if packing nut O-ring (2) is in place inside the valve body groove.
 - d. Do not attempt to repair the valve body if damaged.
2. Components
 - a. Inspect all parts visually for wear, damage. Replace parts as necessary. In case of damage to any part of upper spindle assembly (4)/packing nut assembly (5), replace with new packing nut/upper spindle assembly.
 - b. Inspect lower stem (3) threads and soft seating for any sign of wear / damage. Replace if necessary.
 - c. Inspect PRD (5/9) [if installed] for any damage.
 - d. For parallel inlet connection, replace inlet O-ring (8) once valve is removed from the cylinder.


Assembly of Valve

1. Fit packing nut O-ring (2) inside the groove provided in the valve body (1) just below the packing nut threads.
2. Place the lower stem assembly (3) into the valve body. Engage the upper stem square with the lower stem square and screw in packing nut / upper spindle assembly (4) into the valve body by rotating the upper stem flats. This will also drive the lower stem assembly to rest with the valve body seat.
3. Clamp valve body in bench vice between nylon clamps. Tighten packing nut assembly at 50 Nm in clockwise direction using 13/16" socket wrench.
4. Use 9.50 mm socket wrench or HEX box wrench (for CG-1) or T-30 TORX star tool (for CG-4) for tightening PRD assembly (5/9), if applicable. For CG-1 device, tighten at 16 Nm and for CG-4 device, tighten at 7 Nm in clockwise direction.
5. For parallel inlet connection, fit flange O-ring (6) in the groove provided above the inlet thread, fit the flange ring (7) below the flange O-ring and place the inlet O-ring (8) so that it rests against the flange ring.

NOTE Refer "Material of construction and assembly arrangement" page to identify the part No. given in the bracket.



Series BPTN-12 & BPKN-12



Options – Toggle *
Short toggle (L – 34.5 mm)
Long toggle (L – 44.5 mm)

A - Other inlet, outlet & dip tube connections are available as per customer requirement



List of Approved Gases

Series BPTN-12 & BPKN-12

Sl. No.	Name of Gas ^A	Chemical Formula	Outlet Connection	
			ISO 407	CGA V-1
01	Air	-	FIG 17	CGA 950
02	Carbon dioxide (with & without draw-off) & Carbon dioxide / Oxygen (CO ₂ >7%)	CO ₂	FIG 16	CGA 940
03	Helium & Helium / Oxygen (O ₂ ≤20%)	He	FIG 15	CGA 930
04	Nitrogen	N ₂	FIG 18	CGA 960
05	Nitrous Oxide (with & without draw-off)	N ₂ O	FIG 13	CGA 910
06	Oxygen	O ₂	FIG 9	CGA 870
07	Oxygen / Carbon dioxide (CO ₂ ≤7%)	O ₂ + CO ₂	FIG 10	CGA 880
08	Oxygen / Helium (He ≤80%)	O ₂ + He	FIG 11	CGA 890
09	50% Oxygen / 50% Nitrous Oxide (47.5% <N ₂ O <52.5%)	O ₂ + N ₂ O	FIG 8	CGA 965
10	Non-flammable, non-corrosive gas mixtures labelled as drugs or medical devices and not having another connection assignment	-	FIG 4 (AJ)	CGA 973

Don't let your profit leak through your cylinders



6 Continents
65+ Countries



tekno valves
driven by excellence



Your safety is valued

Tekno Valves

Natun Rasta, Bilkanda, 24 Parganas (N), Kolkata

INDIA

+91 33 25956767

www.teknovalves.com

International Distributors

USA & CANADA

Tekno Valves North America, Inc.

+1 (225) 330 - 6590

www.tvnainc.com

EUROPE

GBP Gas Business Partner GmbH

+49 (0)6468-917 99 52

www.gas-business-partner.com